



# PROTEIN LEVEL OF ENR LATEX

TEST METHOD : ASTM D5712-05

Standard Test Method for Analysis of Aqueous Extractable Protein in Natural Rubber and Its Products Using the Modified Lowry Method

SAMPLE	LOT / BATCH NO.	Average ( $\mu\text{g}$ Protein/g dry rubber)
Natural concentrate latex - HA GRADE	5/2040	1030.65
Natural concentrate latex - LATZ GRADE	21/2049	1082.68
Natural concentrate latex- LOW PROTEIN GRADE	LP 001	468.89
Double Centrifuge natural concentrate latex -DC GRADE	DC 007	237.97
ENR50 LATEX	B390	145.48
ENR25 LATEX	B552	142.04

TESTED BY DEPARTMENT OF BIOCHEMISTRY, FACULTY OF SCIENCE, PRINCE OF SONGKLA UNIVERSITY



**MUANG MAI GUTHRIE**

Public Company Limited.